



## Tree and plant root samples for RNA-Seq

How to dig, clean and preserve root tissue (liquid nitrogen or RNAlater) so the RNA stays intact for RNA-Seq, with quality targets for difficult, polyphenol-rich roots.

Roots contain polysaccharides, polyphenols, soil inhibitors and active RNases, so careful sampling and fast preservation are critical.

### 1. Plan the experiment

- Choose the root type for your question: fine roots (under 2 mm) for nutrient uptake and stress; lateral roots for branching; tap roots for growth; root tips (1–2 cm) for cell division and gene expression.
- At least **3 biological replicates** (independent trees or plants) per treatment, at the same developmental stage and the same time of day.
- Avoid decayed, dark brown or black, or damaged roots.

### 2. Prepare before going to the field

- Liquid nitrogen and cryovials (or RNAlater and RNase-free 2 mL tubes), RNase-free forceps and scissors, gloves, cool box with ice packs, 70% ethanol and an RNase decontaminant (for example RNaseZAP), labels and waterproof marker.
- Label tubes in advance: sample ID, species, location, root type, treatment, date and time, replicate.
- Clean tools between replicates; new gloves for each sample.

### 3. Dig and collect

- Remove soil gently with sterile brushes or forceps; do not pull or crush roots.
- Select healthy roots and cut 2–5 cm segments, or 1–2 cm white, actively growing tips.
- Collect **100–300 mg** fresh root tissue per sample (more for low-yield species).

### 4. Remove soil quickly

Shake off loose soil, rinse briefly in RNase-free water if needed and blot dry. Do not wash for long, as this changes stress-gene expression. Aim for **less than 2–3 minutes** from digging to preservation.

### 5. Preserve immediately

#### Option 1: Liquid nitrogen (best)

Put roots straight into cryovials and snap-freeze. Keep in liquid nitrogen or on dry ice during fieldwork, then store at –80 °C.

#### Option 2: RNAlater (when liquid nitrogen is not available)

1. Cut roots into 0.5–1 cm pieces and put them into RNAlater within about 1 minute.
2. Use 5–10 volumes of RNAlater: about 1 mL per 100 mg tissue. Every piece must be fully covered.
3. Close tightly and invert a few times.
4. Keep at 4 °C for 12–24 hours so RNAlater soaks in, then move to –80 °C.

| Storage in RNAlater | Up to                          |
|---------------------|--------------------------------|
| 25 °C               | 24 hours                       |
| 4 °C                | 1 month                        |
| –20 °C              | Several months                 |
| –80 °C              | More than a year (recommended) |

## 6. Transport

Ship in an insulated box with ice packs (RNAlater samples) or dry ice (frozen samples). Keep below 30 °C and out of the sun. Avoid freeze–thaw.

## 7. Root-specific challenges

| Challenge                 | Effect                           | What to do                                  |
|---------------------------|----------------------------------|---------------------------------------------|
| Polysaccharides           | Low yield                        | CTAB/PVP-based extraction                   |
| Polyphenols               | RNA oxidation                    | β-mercaptoethanol and PVP during extraction |
| Soil                      | Inhibitors, microbial RNA        | Clean roots well before preserving          |
| RNases, slow preservation | Degradation, biased profiles     | Preserve within 1 minute                    |
| Rough handling            | Wound-response genes switched on | Handle gently                               |

## 8. Quality targets for extracted RNA

| Parameter             | Target                           |
|-----------------------|----------------------------------|
| Concentration / total | At least 100 ng/μL and 1 μg      |
| A260/280              | 1.9–2.1                          |
| A260/230              | Above 2.0                        |
| RIN                   | 7 or higher (preferably above 8) |
| DNA, protein, phenol  | Not detectable                   |

Before extraction, pour off excess RNAlater and blot the tissue; carried-over RNAlater can interfere with some kits.

**Questions before sending samples?** WhatsApp / phone +880 1308-222824 · bdgenome@yahoo.com · [bdgenomes.com/guides/plant-root-rna/](http://bdgenomes.com/guides/plant-root-rna/)